


## Filtek™ P60 Posterior Restorative – Fast Cure Time and Excellent Performance



### Key Features

- Superior Packability
- Fast Curing Time
- Low Shrinkage
- Excellent Wear Resistance and Strength
- Long-term Clinical Experience



The Dental Advisor awarded Filtek P60 "Best of 2003 Preferred Products" in the category "Packable".

(Vol. 21, No. 1, January 2004)

With its speed and handling, it's way out front among posterior restoratives. And its stylish beauty makes it a winner with patients too.

### Filtek™ P60 Posterior Restorative offers:

**Superior Packability** – Packs well without sticking to instrument. Excellent marginal adaptation and interproximal contact.

**Easy Sculpting** – With exceptional resistance to slumping.

**Fastest Placement Technique** – Takes only 20 seconds to cure an increment of 2.5mm.


### Low Polymerization

**Shrinkage** – Less stress on the tooth, less sensitivity potential for the patient.

### Excellent Clinical

**Performance** – Filtek P60 was awarded five plus by The Dental Advisor (Vol. 21, No. 4, May 2004) due to the clinical quality of the restorations after five years in clinical service.

Clinical quality from 3M ESPE Filtek P60 over a five years recall period



Restorations were evaluated on a 1-5 rating scale: 1 = poor, 2 = fair, 3 = good, 4 = very good, 5 = excellent

# Broaden Your Experience

**3M ESPE offers a full set of materials and devices for composite restorations – from adhesives to finishing and polishing materials**

The complete programme in an overview: the 3M ESPE Composite Navigator provides orientation and helps to select the best fitting option to your requirements. For further information on products and product families being part of the 3M ESPE Composite Navigator, please contact your 3M ESPE office.

<b>3M ESPE</b>				
Navigation	Composite Navigator for			
	Adhesive Systems	Composite Restoratives	Light Curing Devices	Finishing & Polishing
Overview				
Detailed Information	 	  		
		 		

## Ordering Information

### Item No.

**4720** **Filtek™ P60 Posterior Restorative Introductory Kit**  
 · 3 Syringes Filtek P60, 4 g (A3, B2, C2) each  
 · 1 Bottle Adper Scotchbond 1 XT Adhesive 3ml  
 · 1 Syringe Scotchbond Etchant 3 ml  
 · Accessories, Shade Guide

**4720 A3** **Filtek™ P60 Refill Kits**  
 · 1 Syringe 4g, Shade A3  
**4720 B2** · 1 Syringe 4g, Shade B2  
**4720 C2** · 1 Syringe 4g, Shade C2



## 3M ESPE

3M ESPE AG · ESPE Platz  
 82229 Seefeld · Germany  
 E-mail: [info3mespe@mmm.com](mailto:info3mespe@mmm.com)  
 Internet: <http://www.3mespe.com>

3M, ESPE and Filtek are trademarks of 3M or 3M ESPE AG.