


**coltosol<sup>®</sup> F**

*Temporary Fillingmaterial*




**Economical & reliable**

**COLTENE**

# coltosol® F

## Simple Application, reliable effect



### Easy application

- Creamy consistency for easy adaptation in the cavity
- Good retention in the cavity
- No mixing, no waste of material
- Simple setting by the action of saliva
- Easy removal in large portions
- Long shelf life



### Reliable Security

- Fluoride release
- Radiopacity for the perfect control of the treatment result
- Perfect marginal seal
- Eugenol-free, no effect on the curing of composites
- Ideal surface hardness for provisional fillings in posteriors as well as for root canal sealings

### Instant Application

After rolling the dispensed Coltosol F into a ball with your fingers, simply press it into the prepared cavity.

### Easy removal

Coltosol F can be easily removed in large pieces using a bur or an excavator without leaving any residues.

### The choice is yours

The required quantity of Coltosol F is either dosed with the syringe or taken directly from the glass jar with an instrument. Quick, easy and clean.

### Economical: Optimal value for money



### Order Information:

Coltosol F in syringes:

**5930 Coltosol F**  
5 syringes of 8 g each

Coltosol F in glass jar:

**5906 Coltosol F Single pack**  
1 glass jar of 38 g

**5911 Coltosol F Eco Pack**  
3 jars of 38 g each

www.coltene.com  
© 2014 Coltène/Whaledent AG

#### Coltène/Whaledent AG

Feldwiesenstrasse 20  
9450 Altstätten / Switzerland  
Tel +41 71 757 5300  
Fax +41 71 757 5301  
info.ch@coltene.com

#### Coltène/Whaledent Inc.

235 Ascot Parkway  
Cuyahoga Falls, Ohio 44223 / USA  
Tel +1 330 916 8800  
Fax +1 330 916 7077  
info.us@coltene.com

#### Coltène/Whaledent Ltd.

The President Suite  
Kendal House, Victoria Way  
Burgess Hill, West Sussex  
RH15 9NF / UK  
Tel +44 1444 235486  
Fax +44 1444 870640  
info.uk@coltene.com

 **COLTENE**